

## Certificate of Conformance

Product: Agilent High Sensitivity Protein 250 Ladder  
Included in Part No: 5067-1575, 5067-1578


Part No: G2938-85003  
Lot No: 2307

Note: The reagents in this kit were prepared using ISO 9001:2015 quality processes. **For Research Use Only. Not for use in diagnostic procedures.**

### Kit Information

During verification, the following characteristics were tested for conformance to internal specifications:

Item	Method	Result
<b>High Sensitivity Protein 250 Ladder</b>		
Concentration	Photometrical	Passed <input checked="" type="checkbox"/>
Functional	Electrophoretical	Passed <input checked="" type="checkbox"/>

<b>Kit Expiration Date</b>	13.FEB.2024
Date:	22.02.2023
Approved by:	Jurgita Zilinskiene 
Thermo Fisher Scientific Baltics UAB Vilnius, Lithuania	

