

Certificate of Conformance

July 30, 2013

Justin Suchanec
Industrial Valve & Fitting
Justin_suchanec@ivfinc.com

This is to confirm the following regarding the purchase of the product(s) below:

Supelco Product No: 20409
Product Name: Glass Wool
Lot Number: 40113A
Quantity: 450
Order Number: 3009192452
PO: 0028124

Properties:

packaging

pkg of 10 g

Description:

Pesticide Grade (Silanized) for applications that involve active analytes, such as organochlorine and/or organophosphorous pesticides

Wool plugs are often used in GC inlet liners to improve sample vaporization, and/or to keep non-volatile material from entering the column. They can also be used in packed GC columns, solvent desorption tubes, thermal desorption tubes, and purge traps to retain adsorbent beds.

The material shipped meets our quality control specifications for identity. For further questions, please contact Technical Service by phone or at techserv@sial.com.

Sincerely,

Jeanne

A. Jeanne Harrell/ Technical Service Scientist
Sigma-Aldrich Corporation
3050 Spruce Street/ Saint Louis, MO 63103 / USA
(800) 325-5832, X3701