



CERTIFICATE OF CONFORMANCE

Agilent Product Name: P3311 pH triode combination electrode

Agilent Part No.: P3311 in Kit 5190-3990

Certificate Date: 29 Aug. 2012

Quality Statement

This product is certified and guaranteed to be fully compatible with the Agilent pH meters. All lots of this product are manufactured under an ISO 9001 system. This product was designed and manufactured to a number of critical specifications to ensure the perfect fit with instrumentation.

Quality Test Results

Appearance		Pass
Slope		Pass
Zero		Pass
Response Time		Pass
Temperature Error		Pass

Shelf Life

This product has been determined to have a shelf life of 12 months from the date of manufacture. The product will have best life and performance if used and stored in accordance with the operating guide.

Certified by:

A handwritten signature in black ink that reads "Neil Noronha".

Neil Noronha
Chemistries and Supplies Division
Quality Manager