



Agilent Technologies

CERTIFICATE OF CONFORMANCE

Agilent Product Name: P3213 pH combination electrode

Agilent Part No.: P3213 in kit 5190-3992

Certificate Date: 26 Mar. 2015

Quality Statement

This product is certified and guaranteed to be fully compatible with the Agilent pH meters. All lots of this product are manufactured under an ISO 9001 system. This product was designed and manufactured to a number of critical specifications to ensure the perfect fit with instrumentation.

Quality Test Results

| | | |
|---------------|--|-------------|
| Appearance | | Pass |
| Slope | | Pass |
| Zero | | Pass |
| Response Time | | Pass |

Shelf Life

This product has been determined to have a shelf life of 12 months from the date of manufacture. The product will have best life and performance if used and stored in accordance with the operating guide.

Certified by:

Neil Noronha
Chemistries and Supplies Division
Quality Manager