



**Agilent Technologies**

## CERTIFICATE OF CONFORMANCE

**Agilent Product Name:** P3214 pH combination electrode

**Agilent Part No.:** 5190-3993

**Certificate Date:** 21 Jul. 2014

### Quality Statement

This product is certified and guaranteed to be fully compatible with the Agilent pH meters. All lots of this product are manufactured under an ISO 9001 system. This product was designed and manufactured to a number of critical specifications to ensure the perfect fit with instrumentation.

### Quality Test Results

|               |  |             |
|---------------|--|-------------|
| Appearance    |  | <b>Pass</b> |
| Slope         |  | <b>Pass</b> |
| Zero          |  | <b>Pass</b> |
| Response Time |  | <b>Pass</b> |

### Shelf Life

This product has been determined to have a shelf life of 9 months from the date of manufacture. The product will have best life and performance if used and stored in accordance with the operating guide.

Certified by:

Neil Noronha  
Chemistries and Supplies Division  
Quality Manager