



CERTIFICATE OF CONFORMANCE

Agilent Product Name: D6111 DO probe

Agilent Part No.: D6111 in kit 5190-3997

Certificate Date: 20 Jun. 2012

Quality Statement

This product is certified and guaranteed to be fully compatible with the Agilent pH meters. All lots of this product are manufactured under an ISO 9001 system. This product was designed and manufactured to a number of critical specifications to ensure the perfect fit with instrumentation.

Quality Test Results

Appearance		Pass
Full-Scale Current		Pass
Zero-Oxygen Current		Pass
Response Time		Pass
Temperature Error		Pass

Shelf Life

This product has been determined to have a shelf life of 15 months from the date of manufacture. The product will have best life and performance if used and stored in accordance with the operating guide.

Certified by:

A handwritten signature in black ink that reads "Neil Noronha".

Neil Noronha
Chemistries and Supplies Division
Quality Manager