



## CERTIFICATE OF CONFORMANCE

**Agilent Product Name:** D6111 DO probe

**Agilent Part No.:** D6111 in kit 5190-3997

**Certificate Date:** 28 Jan. 2014

### Quality Statement

This product is certified and guaranteed to be fully compatible with the Agilent pH meters. All lots of this product are manufactured under an ISO 9001 system. This product was designed and manufactured to a number of critical specifications to ensure the perfect fit with instrumentation.

### Quality Test Results

|                     |  |             |
|---------------------|--|-------------|
| Appearance          |  | <b>Pass</b> |
| Full-Scale Current  |  | <b>Pass</b> |
| Zero-Oxygen Current |  | <b>Pass</b> |
| Response Time       |  | <b>Pass</b> |
| Temperature Error   |  | <b>Pass</b> |

### Shelf Life

This product has been determined to have a shelf life of 15 months from the date of manufacture. The product will have best life and performance if used and stored in accordance with the operating guide.

Certified by:

Neil Noronha  
Chemistries and Supplies Division  
Quality Manager