

Agilent CrossLab Product Name: Agilent CrossLab Non-refillable pH Combination Electrode with PC Body

Agilent CrossLab Part No: 8010-0972

Certificate Date: 29 Dec. 2014

Quality Statement

This product is certified and guaranteed to be fully compatible with common pH meters*. All lots of this product are manufactured under an ISO 9001 system. This product was designed and manufactured to a number of critical specifications to ensure the perfect fit with instrumentation.

*For a complete list of compatible pH meters, please see the accompanying compatibility guide.

Quality Test Results

Appearance		Pass
Slope		Pass
Zero		Pass
Response Time		Pass

Shelf Life

This product has been determined to have a shelf life of 12 months from the date of manufacture. The product will have best life and performance if used and stored in accordance with the operating guide.

Certified by:



Neil Noronha
Chemistries and Supplies Division
Quality Manager