

Agilent CrossLab Product Name: Agilent CrossLab Refillable pH Combination Electrode with Glass Body

Agilent CrossLab Part No: 8010-0985

Electrode Model number: 8010-0971

Certificate Date: 26 Mar. 2015

Quality Statement

This product is certified and guaranteed to be fully compatible with common pH meters*. All lots of this product are manufactured under an ISO 9001 system. This product was designed and manufactured to a number of critical specifications to ensure the perfect fit with instrumentation.

*For a complete list of compatible pH meters, please see the accompanying compatibility guide.

Quality Test Results

Electrode Model number: 8010-0971	Appearance	Pass
	Slope	Pass
	Zero	Pass
	Response Time	Pass

Shelf Life

This product has been determined to have a shelf life of 12 months from the date of manufacture. The product will have best life and performance if used and stored in accordance with the operating guide.

Certified by:



Neil Noronha
Chemistries and Supplies Division
Quality Manager