



5067-5601



Declaration of Conformity to Manufacturing Specification

Agilent part number:

5067-5601

Product name:

TapeStation Test Tape

Storage conditions:

Room Temperature

Shipping conditions:

Uncontrolled Temperature

Tape ID:

01-S034-231221-01-000038

Expiration date:

2026-12-20

Release date:

2023-12-21

Approved by:

Sebastian Pallmann

Sebastian Pallmann  
Senior Quality Engineer  
Agilent Technologies  
Deutschland GmbH

Date:

2023-12-21

This product was manufactured and verified in accordance with an ISO 9001 certified quality management system. The manufacture, inspection, and testing were performed by trained personnel. The equipment and materials used in the manufacturing process are identifiable and traceable. All associated test equipment is maintained and calibrated.

The following product characteristics were tested for conformity to manufacturing specifications.

Product characteristic	Method	Result
ScreenTape integrity	Visual inspection	Pass
Barcode integrity	Read test	Pass
Fiducial features	Position check	Pass
Whiteness check	Spectral analysis	Pass
Electric features	Continuity test	Pass

For Research Use Only. Not for use in diagnostic procedures.

Agilent Technologies  
Deutschland GmbH  
Hewlett-Packard-Str. 8  
76337 Waldbronn, Germany

Made in Germany

For global support: [www.agilent.com/en/contact-us/page](http://www.agilent.com/en/contact-us/page)

Template Number: 5067-5601 - DOC MS  
Template Revision: B.00