

Certificate of Conformity to Manufacturing Specifications

FP 165kb Ladder, 35uL

**Agilent Part
Number:**

FP-7002-U035

**Lot
Number:**

0006778192

This analytical reference material was manufactured and verified in accordance with an ISO 9001:2015 certified quality assurance system. The manufacture, inspection, and testing have been performed by trained personnel. The instrumentation and materials used in the qualification process are identified and traceable. All associated test instruments are maintained and calibrated.

These characteristics have been tested on a per lot basis.

Product Characteristic	Method	Result
Appearance	Visual	Pass
Concentration - Fluorometric	Fluorometric	Pass
Purity	Electrophoretical	Pass
Separation Performance	Electrophoretical	Pass

Date of release: Jan 8, 2024

Date of expiration: Jan 1, 2025



Brad Piatt
QA/QC Lead

For Research Use Only. Not for use in diagnostic procedures.

Made in the USA
Agilent Technologies, Inc
2450 SE Oak Tree Ct
Ankeny IA 50021 US

For support within the US or Canada: 1-800-227-9770
For global support: www.agilent.com/en/contact-us/page