

## Certificate of Conformity to Manufacturing Specifications

Intercalating Dye, 30uL

**Agilent Part  
Number:**

DNF-600-U030

**Lot  
Number:**

0006777880

This analytical reference material was manufactured and verified in accordance with an ISO 9001:2015 certified quality assurance system. The manufacture, inspection, and testing have been performed by trained personnel. The instrumentation and materials used in the qualification process are identified and traceable. All associated test instruments are maintained and calibrated.

These characteristics have been tested on a per lot basis.

Product Characteristic	Method	Result
Appearance	Visual	Pass
Functional Performance	Electrophoretical	Pass

**Date of release:** Jan 9, 2024

**Date of expiration:** Dec 26, 2024



---

Brad Piatt  
QA/QC Lead

For Research Use Only. Not for use in diagnostic procedures.

**Made in the USA**  
Agilent Technologies, Inc  
2450 SE Oak Tree Ct  
Ankeny IA 50021 US

**For support within the US or Canada: 1-800-227-9770**  
**For global support: [www.agilent.com/en/contact-us/page](http://www.agilent.com/en/contact-us/page)**