

## Certificate of Conformity to Manufacturing Specifications

**NGS Fragment Separation Gel, 240mL**

**Agilent Part  
Number:**

DNF-240-0240

**Lot  
Number:**

0006777095

This analytical reference material was manufactured and verified in accordance with an ISO 9001:2015 certified quality assurance system. The manufacture, inspection, and testing have been performed by trained personnel. The instrumentation and materials used in the qualification process are identified and traceable. All associated test instruments are maintained and calibrated.

These characteristics have been tested on a per lot basis.

Product Characteristic	Method	Result
Appearance	Visual	Pass
Conductivity	Electrochemical	Pass
Microbial Testing	Enumeration	Pass
Separation Performance	Electrophoretical	Pass
Viscosity	Viscometry	Pass
pH	Electrochemical	Pass

For Research Use Only. Not for use in diagnostic procedures.

**Made in the USA**

Agilent Technologies, Inc  
2450 SE Oak Tree Ct  
Ankeny IA 50021 US

**For support within the US or Canada: 1-800-227-9770**  
**For global support: [www.agilent.com/en/contact-us/page](http://www.agilent.com/en/contact-us/page)**

**Date of release:** Jan 10, 2024

**Date of expiration:** Dec 18, 2024



---

Brad Piatt  
QA/QC Lead

For Research Use Only. Not for use in diagnostic procedures.

**Made in the USA**  
Agilent Technologies, Inc  
2450 SE Oak Tree Ct  
Ankeny IA 50021 US

**For support within the US or Canada: 1-800-227-9770**  
**For global support: [www.agilent.com/en/contact-us/page](http://www.agilent.com/en/contact-us/page)**