

## Declaration of Conformity to Manufacturing Specifications

Product: Agilent High Sensitivity DNA Chips

Part No: G2938-68000

Included in Part No: 5067-4626

Lot No: BZ18BK50

Note: The chips were prepared using ISO 9001 quality processes. **For Research Use Only. Not for use in diagnostic procedures.**

### Kit Information

During specification verification, the following characteristics were tested on a statistically significant number of chips:

Specification	Method	Result
Channel Depth	Profilometer	Passed <input checked="" type="checkbox"/>
Optical Inspections	Visual Inspections	Passed <input checked="" type="checkbox"/>

Kit Expiration Date:	18.DEC.2024
Date: 10.JAN.2024	
Signed: <b>S.Knapp</b>	
PerkinElmer, Inc. Hopkinton, MA 01748	



**Agilent Technologies**